

OPU ProbeHolder

User Instruction Manual



OPU ProbeHolder

After sales correspondence

All correspondence relating to the after sales service of Podoblock B.V.

to be addressed to:

NL Customers

Podoblock B.V.

Verlengde Nijverheidsweg 10

9482 WT Tynaarlo

The Netherlands

Tel: +31 50 8200 257

Fax: +31 50 8200 212

Overseas Customers

Contact your local distributor. In case of doubt please contact

Podoblock B.V.

For further information visit www.podoblock.com

Trade marks

Podoblock™ is a trademark of Podoblock B.V.

Warning!

Please read this preliminary information carefully and note ALL of the warnings, cautions and safety notes contained within these instructions for use before using the OPU ProbeHolder. Keep these Instructions For Use (IFU) close at hand at all times for reference.

The ProbeHolder must only be used as detailed within this “handbook”. Failure to do so could result in injury to patients or users, or (also) damage to the equipment. Podoblock does not take any responsibility for improper handling of the OPU ProbeHolder.

- The instruction for use should always be readily accessible for reference and be referred to for details of the OPU ProbeHolder.
- Sections marked with the symbol provide guidance to possible actions that could lead to injury of the patient or user, or damage to the ProbeHolder.
- Ensure all warnings in the handbook are strictly complied with.
- Section 3 Operation provides information detailing how to use the ProbeHolder.
- The OPU ProbeHolder is compatible with most Ultrasound source. Never use the OPU ProbeHolder without first carefully reading this manual. Correct handling may prevent damage to the product or harm the user. Misuse, abuse of the OPU ProbeHolder will not be covered by the given warranty, neither will warranty cover the damage because of accident or high impact.

OPU Probe Holder for equine ultrasounds and is compatible with most ultrasound device brands on the market.

With the ProbeHolder, you have an innovation that makes equine OPU ultrasound processes simpler, faster and more comfortable for both vet and horse.

The unique construction, material and size of this innovative product are the attributes that make it valuable. It comes equipped with an integrated needle chamber which can be (optional) supplied with a stainless steel needle guide.

Placing your Ultrasound probe is quick and easy and requires only one locking screw through the handle with integrated magnets secure it all.



Packaging

The OPU ProbeHolder is packaged and transported in a hardcover protective case and cushioned with protective foam.



3


Operating Instructions

1. Check carefully that all parts are included in the hardcase:

- ProbeHolder
- Needle Guide (Optional)

Follow the instructions below

2. Unscrew the red knob screw by turning it counter-clockwise



3. Slide the complete handle grip backwards




4. Open the ProbeHolder by placing finger in the desired notches. Both halves are secured by 2 magnets



5. Carefully insert the needle through the needle guide and close the probe holder by clapping the two halves together and screw back on the red screw to lock the probe holder securely



6. The black removable handgrip can be placed in three (x3) positions (removed by turning round by hand)





Podoblock B.V.

Verlengde Nijverheidsweg 10

9482 WT Tynaarlo

The Netherlands

info@podoblock.com

0031 50 8200 257

fax: 0031 50 8200 212